

# Vira-Care

a Broad Spectrum Disinfectant-Cleaner

*Takes care of all your cleaning needs!*

## AREAS FOR USAGE

## DILUTION RATIOS

(preferably with warm water)

### Isolation/Quarantine / Problem Areas / Surgical Tools / Equipment

Remove all solids prior to soaking or washing down – use on floors, walls and all hard surfaces.

Allow plenty of coverage and sufficient contact time. Leave to dry.

**1 : 5**

Soak surgical tools and equipment for approx 10 – 15 mins - **DO NOT RINSE - ALLOW TO DRY**

### ENCLOSURES & FLOORS (both Inside and Outside)

Floors (heavily soiled)

**1 : 50**

We recommend all solids are removed prior to washing and the product is diluted between:

**1 : 100**

ALLOW TO DRY.

Floors (soiled)

**1 : 100**

We recommend all solids are removed prior to washing and the product is diluted between:

**1 : 200**

ALLOW TO DRY.

Changeovers and Weekly Washes dilute between:

**1 : 5**

and always ALLOW TO DRY.

**1 : 100**

General Cleaning – Floors, walls, surfaces etc

**1 : 200**

### HARD SURFACES (Ideal for all hard surfaces inc wood, metal, glass, plastic, rubber and polymer)

General, including worktops, cupboards, sinks, shelves, holding areas, wooden decking, paths, glass

**1 : 50**

surfaces, Trays and Crates, Coops and other animal housing types. Dependant on the level of

**1 : 200**

soiling/deposits (which must be removed first) dilute between -

### BOWLS / CONTAINERS / TOOLS / UTENSILS

We recommend you wash first to remove solids. Soak for approx 10 - 15 minutes, scrub to remove

**1 : 50**

any stubborn deposits, if necessary, then rinse well and allow to dry.

**1 : 100**

### Footbaths:

Ideal for footbaths. Make up a fresh batch regularly, particularly when water becomes muddied.

**1 : 50**

### Multi-Purpose & Effective against: Bacteria, Viruses, Fungi, Algae, Yeast, Protozoa and Mould

Container Size to be made-up	Dilution Rqud	Pump Depressions	ml	Container Size to be made-up	Dilution Rqud	Pump Depressions	ml
5 litre	1:50	4	100ml	10 litre	1:100	4	100ml
	1:100	2	50ml		1:200	2	50ml
	1:200	1	25ml				
15 litre	1:100	6	150ml	20 litre	1:100	8	200ml
	1:200	3	75ml		1:200	4	100ml
25 litre	1:100	10	250ml	30 litre	1:100	12	300ml
	1:200	5	125ml		1:200	6	150ml

Each pump depression should provide approximately 25ml of product.

**Vets, Kennels, Catteries, Stables, Livestock and Avian Environments**



# GHS direct

For further information or technical support please contact:  
01268 513200 [sales@ghs-direct.com](mailto:sales@ghs-direct.com) [www.ghs-direct.com](http://www.ghs-direct.com)

